

Backup Masks for Medical Personnel

Currently in the situation that we are in some people are now making masks for medical personnel that meet requirements for the CDC. Below are the requirements for the mask, materials needed, and the procedures.

Requirements:

- Must be tight woven fabric
- Must have a pocket for a N95 filter
- Must have a contact person at facility to organize a drop time and location
- Must be sanitized at location of use.

Materials: (tutorial to be posted soon)

- 9" x 15" piece of 100% quilting cotton (solid or print)
- (2) 7.5 "pieces of 1/4" or 5/8" braided plastic
- 6" pipe cleaner
- Sewing Machine with matching thread



- Iron/
Iron-
ing



- Scis-

sors/

Ruler

- Straight Pins or Wonderclips

Tutorial is below:

https://vimeo.com/user109783413/download/399451183/4de0f322d1?fbclid=IwAR3ZCyxa6h9qq8YiV-CANUKhO_q-wmIsePE4Zs__WKQXH99m0121jywGHgjo

Procedures:

- Gain approval by hospital board for use of these masks (we have it)
- Arrange a drop of location where people can still stay a safe distance away without overloading hospital
- Have a clear bin with cover for masks
- Once bin is brought to facility an appointed hospital employee will retrieve the bin with a lid and will have it brought to laundry to clean.